

## 3.A3.E Worksheet E – Nested Amplification of Nucleic Acid by PCR for the Confirmation of *Piscirickettsia salmonis*

Case Number \_\_\_\_\_

Date \_\_\_\_\_

PCR Reagent	Lot#	Final Concentration	Stock Concentration**	Volume per Reaction (µL) (to total 50 µL)	Volume for _____ samples
d-H <sub>2</sub> O*		Add to total 48 µL		38.0	
10XBuffer		1X	10X	5.0	
MgCL <sub>2</sub>		1.5 mM	50 mM	1.5	
dNTP's		0.2 mM	10 mM	1	
(+)Primer		1 mM	1 mM	1	
(-)Primer		1 mM	1 mM	1	
TAQ		2.5 units/Rx	5U/µL	0.5	
DNA <sup>±</sup> (round 1)		-	-	2 µL	-

\*Add water to Master Mix first, TAQ last. \*\*Change "Stock Concentration" parameters as necessary. Different reagent sources supply varying stock concentrations. †Do not add DNA template until Master Mix reaction tubes have been removed from the reagent mixing (MM) area.

### Primer Sets for *P. salmonis* 2<sup>nd</sup> (Nested) Round Amplification

PS2S (223F) (round 2 forward)	5' - CTA-GGA-GAT-GAG- CCC-GCG-TTG -3'
PS2AS (690R) (round 2 reverse)	5' - GCT-ACA-CCT-GAA-ATT-CCA-CTT -3'

Gel Concentration	Weight of agarose (grams)	Volume of Buffer (mL)

### Gel Template (Sample Placement Map)

Ladder Brand / Lot # \_\_\_\_\_

Loading Buffer Brand / Lot # \_\_\_\_\_

Enter sample ID below corresponding well number:

PCR Products (Loaded LEFT to RIGHT)														
<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>	<u>13</u>	<u>14</u>	<u>15</u>
<u>16</u>	<u>17</u>	<u>18</u>	<u>19</u>	<u>20</u>	<u>21</u>	<u>22</u>	<u>23</u>	<u>24</u>	<u>25</u>	<u>26</u>	<u>27</u>	<u>28</u>	<u>29</u>	<u>30</u>