

## 3.A3.F Worksheet F – Direct Amplification of Nucleic Acid by PCR for the Confirmation of *Piscirickettsia salmonis*

Case Number \_\_\_\_\_

Date \_\_\_\_\_

PCR Reagent	Lot#	Final Concentration	Stock Concentration**	Volume per Reaction (µL) (to total 50 µL)	Volume for _____ samples
d-H <sub>2</sub> O*		Add to total 45 µL		35.0	
10XBuffer		1X	10X	5.0	
MgCL <sub>2</sub>		1.5 mM	50 mM	1.5	
dNTP's		0.2 mM	10 mM	1	
(+)Primer		1 mM	1 mM	1	
(-)Primer		1 mM	1 mM	1	
TAQ		2.5 units/Rx	5U/µL	0.5	
DNA <sup>±</sup>		-	-	5 µL	-

\*Add water to Master Mix first, TAQ last. \*\*Change “Stock Concentration” parameters as necessary. Different reagent sources supply varying stock concentrations. <sup>±</sup>Do not add DNA template until Master Mix reaction tubes have been removed from the reagent mixing (MM) area.

### Primer Sets for *P. salmonis* 2<sup>nd</sup> Direct Amplification

<b>PS2S (223F) (round 2 forward)</b>	5' - CTA-GGA-GAT-GAG- CCC-GCG-TTG -3'
<b>PS2AS (690R) (round 2 reverse)</b>	5' - GCT-ACA-CCT-GAA-ATT-CCA-CTT -3'

Gel Concentration	Weight of agarose (grams)	Volume of Buffer (mL)

### Gel Template (Sample Placement Map)

Ladder Brand / Lot # \_\_\_\_\_

Loading Buffer Brand / Lot # \_\_\_\_\_

Enter sample ID below corresponding well number:

PCR Products (Loaded LEFT to RIGHT)														
<u>1</u>	<u>2</u>	<u>3</u>	<u>4</u>	<u>5</u>	<u>6</u>	<u>7</u>	<u>8</u>	<u>9</u>	<u>10</u>	<u>11</u>	<u>12</u>	<u>13</u>	<u>14</u>	<u>15</u>
<u>16</u>	<u>17</u>	<u>18</u>	<u>19</u>	<u>20</u>	<u>21</u>	<u>22</u>	<u>23</u>	<u>24</u>	<u>25</u>	<u>26</u>	<u>27</u>	<u>28</u>	<u>29</u>	<u>30</u>